



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human ADRB1 Antibody (6B9), FITC

#### OVERVIEW

Reference	ARO-A10776
Clonality	Monoclonal Antibody
Isotype	VHH-hFc

#### PRODUCT INFORMATION

Host Species	Human
Target	Beta-1 adrenoceptor, ADRB1, Beta-1 adrenergic receptor, ADRB1R, B1AR, Beta-1 adrenoceptor
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Beta-1 adrenoceptor, ADRB1, Beta-1 adrenergic receptor, ADRB1R, B1AR, Beta-1 adrenoceptor