

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human CD98/SLC3A2 Antibody (Ab110), FITC
<b>PRODUCT REFERENCE</b>	ARO-A10795
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1, kappa
<b>TARGET</b>	SLC3A2, 4F2hc, Lymphocyte activation antigen 4F2 large subunit, CD98, 4F2 cell-surface antigen heavy chain, Solute carrier family 3 member 2, MDU1, 4F2 heavy chain antigen

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Mus musculus (Mouse)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 0.2% BSA, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-SLC3A2, 4F2hc, Lymphocyte activation antigen 4F2 large subunit, CD98, 4F2 cell-surface antigen heavy chain, Solute carrier family 3 member 2, MDU1, 4F2 heavy chain antigen