



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Rat CD2/LFA-2/OX-34 Antibody (OX-34), FITC

#### OVERVIEW

Reference	ARO-A10832
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	T-cell surface antigen CD2, LFA-2, LFA-3 receptor, OX-34 antigen, T-cell surface antigen T11/Leu-5, CD2, Cd2
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-T-cell surface antigen CD2, LFA-2, LFA-3 receptor, OX-34 antigen, T-cell surface antigen T11/Leu-5, CD2, Cd2