



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human FGF1/aFGF Antibody (SAA0436), FITC

#### OVERVIEW

Reference	ARO-A10845
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	Endothelial cell growth factor, HBGF-1, Fibroblast growth factor 1, FGF1, ECGF, Acidic fibroblast growth factor, aFGF, Heparin-binding growth factor 1, FGF-1, FGFA
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Endothelial cell growth factor, HBGF-1, Fibroblast growth factor 1, FGF1, ECGF, Acidic fibroblast growth factor, aFGF, Heparin-binding growth factor 1, FGF-1, FGFA