



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD340/ERBB2/HER2 Antibody (MAB 2C4), FITC

OVERVIEW

| | |
|-----------|---------------------|
| Reference | ARO-A10864 |
| Clonality | Monoclonal Antibody |
| Isotype | IgG2a, kappa |

PRODUCT INFORMATION

| | |
|---------------------|---|
| Host Species | Mouse |
| Target | p185erbB2, HER2, NGL, Tyrosine kinase-type cell surface receptor HER2, NEU, MLN 19, Proto-oncogene Neu, MLN19, ERBB2, Proto-oncogene c-ErbB-2, CD340, Receptor tyrosine-protein kinase erbB-2, Metastatic lymph node gene 19 protein |
| Form | Liquid |
| Buffer | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300. |
| Stability & Storage | 4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| Aliases /Synonyms | Anti-p185erbB2, HER2, NGL, Tyrosine kinase-type cell surface receptor HER2, NEU, MLN 19, Proto-oncogene Neu, MLN19, ERBB2, Proto-oncogene c-ErbB-2, CD340, Receptor tyrosine-protein kinase erbB-2, Metastatic lymph node gene 19 protein |