



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Normal Monkey (Rhesus) IgG Isotype Control, FITC

#### OVERVIEW

---

Reference	ARO-A10900
Clonality	Polyclonal Antibody
Isotype	IgG

#### PRODUCT INFORMATION

---

Host Species	Rhesus Monkey
Target	IgG, Immune globulin G
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-IgG, Immune globulin G