

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human IL1B/IL1F2 Antibody (SAA0369), FITC
PRODUCT REFERENCE	ARO-A10923
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Interleukin-1 beta, IL1B, Catabolin, IL1F2, IL-1 beta,IL-1 $\hat{1}^2$

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 0.2% BSA, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Interleukin-1 beta, IL1B, Catabolin, IL1F2, IL-1 beta,IL-1 $\hat{1}^2$