



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CXCL11/I-TAC Antibody (SAA0466), FITC

#### OVERVIEW

Reference	ARO-A10991
Clonality	Monoclonal Antibody
Isotype	VHH, mFc

#### PRODUCT INFORMATION

Host Species	Mouse
Target	IP-9, I-TAC, C-X-C motif chemokine 11, ITAC, Small-inducible cytokine B11, SCYB11, Beta-R1, Interferon gamma-inducible protein 9, H174, SCYB9B, CXCL11, Interferon-inducible T-cell alpha chemoattractant
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-IP-9, I-TAC, C-X-C motif chemokine 11, ITAC, Small-inducible cytokine B11, SCYB11, Beta-R1, Interferon gamma-inducible protein 9, H174, SCYB9B, CXCL11, Interferon-inducible T-cell alpha chemoattractant