



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-European white birch Bet v 1/Allergen Bet v I-A Polyclonal Antibody

OVERVIEW

Reference	ARO-A11002
Clonality	Polyclonal Antibody
Isotype	

PRODUCT INFORMATION

Host Species	Rabbit
Target	
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Major pollen allergen Bet v 1-A, Allergen Bet v I-A, Bet v 1-A, BETVIA, BETVI, Betula pendula (European white birch) (Betula verrucosa)