



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-MARV VP40/Marburg VP40 Polyclonal Antibody

#### OVERVIEW

Reference	ARO-A11006
Clonality	Polyclonal Antibody
Isotype	

#### PRODUCT INFORMATION

Host Species	Rabbit
Target	
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Matrix protein VP40, Marburg VP40, mVP40, Membrane-associated protein VP40, VP40, Lake Victoria marburgvirus (strain Musoke-80) (MARV) (Marburg virus (strain Kenya/Musoke/1980))