

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse ST6GALNAC2/SIAT7-B Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11019
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Mouse ST6GALNAC2/SIAT7-B (Thr97-Arg373).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 2 (EC 2.4.3.3) (Gal-beta-1,3-GalNAc alpha-2,6-sialyltransferase) (GalNAc alpha-2,6-sialyltransferase II) (ST6GalNAc II) (ST6GalNAcII) (Sialyltransferase 7B) (SIAT7-B)