

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|----------------------------------------------------------------|
| PRODUCT NAME | Anti-Mouse VSIR/VISTA Polyclonal Antibody |
| PRODUCT REFERENCE | ARO-A11046 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | |
| TARGET | E. coli - derived recombinant Mouse VSIR/VISTA (Ala32-Asn181). |

PRODUCT INFORMATION

| | |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ORIGIN SPECIES | <i>Oryctolagus cuniculus (Rabbit)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-V-type immunoglobulin domain-containing suppressor of T-cell activation (Platelet receptor Gi24) (V-set domain-containing immunoregulatory receptor) (V-set immunoregulatory receptor) |