

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse CD200R1L Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11061
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Mouse CD200R1L (Ile25-Thr220).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Cell surface glycoprotein CD200 receptor 2 (CD200 cell surface glycoprotein receptor-like 2) (CD200 receptor-like 2) (CD200 cell surface glycoprotein receptor-like c) (CD200RLc) (Cell surface glycoprotein OX2 receptor 2)