

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human GFM1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11137
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human GFM1 (Arg526-Asn751).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Elongation factor G1, GFM, GFM1, EFG1, mEF-G 1, Elongation factor G, mitochondrial, EF-Gmt, EFG, hEFG1, Elongation factor G 1, mitochondrial