

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ARID3A Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11160
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human ARID3A (Gly221-Arg351).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Bright, DRIL3, E2F-binding protein 1, ARID domain-containing protein 3A, ARID3A, DRX, E2FBP1, DRIL1, AT-rich interactive domain-containing protein 3A, B-cell regulator of IgH transcription, Dead ringer-like protein 1