

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human INTS13 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11167
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human INTS13 (Ser50-His254).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Protein asunder homolog, C12orf11, Cell cycle regulator Mat89Bb homolog, Integrator complex subunit 13, INTS13, Germ cell tumor 1, Sarcoma antigen NY-SAR-95, GCT1, ASUN