

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human SFTPA1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11208
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human SFTPA1 (Ser130-Phe248).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Collectin-4, Pulmonary surfactant-associated protein A1, SFTP1, SP-A1, COLEC4, SP-A, PSP-A, SFTPA1, SFTPA1B, Alveolar proteinosis protein, SFTPA, PSPA, PSAP, 35 kDa pulmonary surfactant-associated protein