

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human AP1M1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11219
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human AP1M1 (Met1-Pro258).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Clathrin coat-associated protein AP47, Adaptor-related protein complex 1 subunit mu-1, Mu-adaptin 1, Adaptor protein complex AP-1 subunit mu-1, Clathrin coat assembly protein AP47, Golgi adaptor HA1/AP1 adaptin mu-1 subunit, Mu1A-adaptin, AP-1 comple