

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human HMG20B Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11234
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human HMG20B (Tyr77-Gly171).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-HMGX2, BRAF35, HMG domain-containing protein HMGX2, HMG domain-containing protein 2, HMGXB2, Sox-like transcriptional factor, SMARCE1-related protein, HMG box-containing protein 20B, BRCA2-associated factor 35, HMG20B, Structural DNA-binding protein