

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human HIF2A/EPAS1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11237
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human HIF2A/EPAS1 (Leu239-Asn350).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Basic-helix-loop-helix-PAS protein MOP2, bHLHe73, HIF-2-alpha, Member of PAS protein 2, PAS domain-containing protein 2, HIF2A, HLF, EPAS1, Endothelial PAS domain-containing protein 1, EPAS-1, HIF-1-alpha-like factor, Hypoxia-inducible factor 2-alpha