

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human SUZ12 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11252
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human SUZ12 (Lys587-Gly691).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-SUZ12, Polycomb protein SUZ12, Joined to JAZF1 protein, CHET9, KIAA0160, ChET 9 protein, Suppressor of zeste 12 protein homolog, Chromatin precipitated E2F target 9 protein, JJAZ1