

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human EFNA3 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11293
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human EFNA3 (Gly30-Val165).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-EPLG3, EPH-related receptor tyrosine kinase ligand 3, EFL2, LERK-3, EFNA3, Ephrin-A3, LERK3, EFL-2, EHK1-L, EHK1 ligand