

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human LMOD1 Polyclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A11313
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	
<b>TARGET</b>	E. coli - derived recombinant Human LMOD1 (Ser316-Glu497).

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Oryctolagus cuniculus (Rabbit)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-SM-Lmod, Leiomodrin, muscle form, Leiomodrin-1, Smooth muscle leiomodrin, LMOD1, Thyroid-associated ophthalmopathy autoantigen, 64 kDa autoantigen 1D3, 64 kDa autoantigen 1D, 64 kDa autoantigen D1