

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human LAMC1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11335
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human LAMC1 (Gly70-Lys287).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Laminin-7 subunit gamma, Laminin-1 subunit gamma, Laminin subunit gamma-1, Laminin-10 subunit gamma, Laminin-2 subunit gamma, S-laminin subunit gamma, LAMC1, Laminin-11 subunit gamma, Laminin-3 subunit gamma, S-LAM gamma, Laminin-6 subunit gamma, Lam