

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Salmonella typhimurium murD Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11392
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Salmonella typhimurium murD (Pro32-Gly438).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-UDP-N-acetylmuramoylalanine--D-glutamate ligase, 6.3.2.9, D-glutamic acid-adding enzyme, UDP-N-acetylmuramoyl-L-alanyl-D-glutamate synthetase, murD, SL1344_0126, G1W50_13885, G1W53_14860, G1W54_15240, G1W55_14350, G1W56_15745, G1W57_13765, G1W63_1555