

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human KRT81 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11454
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human KRT81 (Ser266-Gly420).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Keratin, hair, basic, 1, KRT81, ghHKb1, Keratin, type II cuticular Hb1, Type-II keratin Kb21, Keratin-81, Type II hair keratin Hb1, MLN 137, K81, ghHb1, Hair keratin K2.9, MLN137, Metastatic lymph node 137 gene protein, KRTHB1