

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|-----------------------------------------------------------|
| PRODUCT NAME | Anti-Human BCL11A Polyclonal Antibody |
| PRODUCT REFERENCE | ARO-A11455 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | |
| TARGET | E. coli - derived recombinant Human BCL11A (Met1-Arg210). |

PRODUCT INFORMATION

| | |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ORIGIN SPECIES | <i>Oryctolagus cuniculus (Rabbit)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-B-cell CLL/lymphoma 11A, EVI-9, CTIP1, EVI9, Zinc finger protein 856, KIAA1809, B-cell lymphoma/leukemia 11A, COUP-TF-interacting protein 1, Ecotropic viral integration site 9 protein homolog, ZNF856, BCL11A, BCL-11A |