

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ACSL1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11462
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human ACSL1 (Gly457-Val698).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-FACL1, LACS2, Long-chain acyl-CoA synthetase 1, Long-chain acyl-CoA synthetase 2, Long-chain-fatty-acid--CoA ligase 1, ACS1, LACS 1, Palmitoyl-CoA ligase 2, Phytanate--CoA ligase, FACL2, Arachidonate--CoA ligase, Long-chain fatty acid-CoA ligase 2, A