

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Anti-Human NCF1 Polyclonal Antibody |
| PRODUCT REFERENCE | ARO-A11486 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | |
| TARGET | E. coli - derived recombinant Human NCF1 (Met1-Val390). |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Oryctolagus cuniculus (Rabbit)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-SH3PXD1A, NCF1, p47-phox, 47 kDa autosomal chronic granulomatous disease protein, Nox-organizing protein 2, NOXO2, SH3 and PX domain-containing protein 1A, NCF-1, 47 kDa neutrophil oxidase factor, NCF-47K, Nox organizer 2, Neutrophil NADPH oxidase fa |