

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human GLI3 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11517
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human GLI3 (Met111-Ser343).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-GLI3 form of 190 kDa, GLI3FL, GLI3-83, Transcriptional activator GLI3, GLI3, GLI3-190, GLI3 C-terminally truncated form, GLI3 full-length protein, GLI3 form of 83 kDa