

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human AP2M1 Polyclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A11575
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	
<b>TARGET</b>	E. coli - derived recombinant Human AP2M1 (Arg170-Cys435).

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Oryctolagus cuniculus (Rabbit)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-Adaptor-related protein complex 2 subunit mu, AP-2 mu chain, CLAPM1, HA2 50 kDa subunit, Plasma membrane adaptor AP-2 50 kDa protein, Adaptin-mu2, AP2M1, KIAA0109, Clathrin coat assembly protein AP50, Adaptor protein complex AP-2 subunit mu, Clathrin