

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human HSD17B1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11604
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human HSD17B1 (Thr4-Pro289).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-EDH17B1, Short chain dehydrogenase/reductase family 28C member 1, EDHB17, Placental 17-beta-hydroxysteroid dehydrogenase, SDR28C1, 20 alpha-hydroxysteroid dehydrogenase, 20-alpha-HSD, HSD17B1, 17-beta-HSD 1, E2DH, Estradiol 17-beta-dehydrogenase 1, E