

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ALDH7A1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11663
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human ALDH7A1 (Met1-Gln539).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ALDH7A1, Delta1-piperideine-6-carboxylate dehydrogenase, Antiquitin-1, P6c dehydrogenase, Alpha-AASA dehydrogenase, Alpha-aminoadipic semialdehyde dehydrogenase, Aldehyde dehydrogenase family 7 member A1, Betaine aldehyde dehydrogenase, ATQ1