

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CCL4L1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11708
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human CCL4L1 (Ala24-Asn92).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-MIP-1-beta, CCL4L, SCYA4L2, Monocyte adherence-induced protein 5-alpha, CCL4L1, Lymphocyte activation gene 1 protein, LAG-1, SCYA4L1, Macrophage inflammatory protein 1-beta, Small-inducible cytokine A4-like, LAG1, C-C motif chemokine 4-like