

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human MYL9 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11732
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human MYL9 (Arg5-Leu163).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Myosin regulatory light chain MRLC1, 20 kDa myosin light chain, Myosin regulatory light chain 2, smooth muscle isoform, MYRL2, MYL9, MLC-2C, LC20, Myosin regulatory light chain 9, MLC2, MRLC1, Myosin RLC, Myosin regulatory light polypeptide 9