

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human POTE ^G Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11756
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human POTE ^G (Leu177-Asp330).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-POTE ^G , POTE14, POTE-14, POTE ankyrin domain family member G, ANKRD26-like family C member 2, Prostate, ovary, testis-expressed protein on chromosome 14, A26C2