

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CAPN3 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11768
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human CAPN3 (Ile602-Ala821).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-CANPL3, New calpain 1, nCL-1, Calpain p94, Calpain L3, CANP3, CAPN3, Calcium-activated neutral proteinase 3, Muscle-specific calcium-activated neutral protease 3, Calpain-3, NCL1, CANP 3