

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human NABP1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11845
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human NABP1 (Asn22-Gly126).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Nucleic acid-binding protein 1, SOSS complex subunit B2, OBFC2A, Single-stranded DNA-binding protein 2, Sensor of ssDNA subunit B2, hSSB2, SSB2, SOSS-B2, NABP1, Oligonucleotide/oligosaccharide-binding fold-containing protein 2A, Sensor of single-stra