

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human APBB1IP Polyclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A11858
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	
<b>TARGET</b>	E. coli - derived recombinant Human APBB1IP (Lys177-Gly420).

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Oryctolagus cuniculus (Rabbit)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-Proline-rich protein 73, RIAM, RARP-1, APBB1-interacting protein 1, Rap1-GTP-interacting adapter molecule, PREL-1, Proline-rich EVH1 ligand 1, Amyloid beta A4 precursor protein-binding family B member 1-interacting protein, PREL1, RARP1, APBB1IP, Ret