

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human GADD45GIP1 Polyclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A11861
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	
<b>TARGET</b>	E. coli - derived recombinant Human GADD45GIP1 (Thr47-Ser222).

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Oryctolagus cuniculus (Rabbit)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-CRIF1, GADD45GIP1, 39S ribosomal protein L59, mitochondrial, PLINP1, PLINP, MRP-L59, CKII beta-associating protein, Papillomavirus L2-interacting nuclear protein 1, PLINP-1, PRG6, CR6-interacting factor 1, MRPL59, p53-responsive gene 6 protein, Growt