

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human UBE2D3 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11870
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human UBE2D3 (Met1-Met147).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Ubiquitin-conjugating enzyme E2-17 kDa 3, (E3-independent) E2 ubiquitin-conjugating enzyme D3, UBC5C, UBE2D3, Ubiquitin-conjugating enzyme E2(17)KB 3, Ubiquitin carrier protein D3, Ubiquitin-protein ligase D3, UBCH5C, Ubiquitin-conjugating enzyme E2