

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human TRIM69 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11874
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human TRIM69 (Thr331-Gln500).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Tripartite motif-containing protein 69, RFP-like domain-containing protein trimless, RNF36, E3 ubiquitin-protein ligase TRIM69, RING finger protein 36, RING-type E3 ubiquitin transferase TRIM69, TRIM69