

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human HECTD1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11884
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human HECTD1 (Lys659-Thr801).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-HECT-type E3 ubiquitin transferase HECTD1, E3 ubiquitin-protein ligase HECTD1, HECTD1, KIAA1131, EULIR, HECT domain-containing protein 1, E3 ligase for inhibin receptor