

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human DPF2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A11894
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human DPF2 (Glu290-Ser391).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Protein requiem, BAF45D, Zinc finger protein ubi-d4, BRG1-associated factor 45D, REQ, DPF2, UBID4, Apoptosis response zinc finger protein, D4, zinc and double PHD fingers family 2