



# ProteoGenix

**ProteoGenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

**ProteoGenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

products@proteogenix.fr

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |                                                         |
|--------------------------|---------------------------------------------------------|
| <b>PRODUCT NAME</b>      | Anti-Human ARL2 Polyclonal Antibody                     |
| <b>PRODUCT REFERENCE</b> | ARO-A12022                                              |
| <b>CLONALITY</b>         | Polyclonal Antibody                                     |
| <b>ISOTYPE</b>           |                                                         |
| <b>TARGET</b>            | E. coli - derived recombinant Human ARL2 (Met1-Asp184). |

### PRODUCT INFORMATION

---

|                                |                                                                                                                                                                                                                                                                                                                     |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ORIGIN SPECIES</b>          | <i>Oryctolagus cuniculus (Rabbit)</i>                                                                                                                                                                                                                                                                               |
| <b>FORM</b>                    | Liquid                                                                                                                                                                                                                                                                                                              |
| <b>BUFFER</b>                  | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300                                                                                                                                                                                                                                                                 |
| <b>STABILITY &amp; STORAGE</b> | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| <b>ALIAS</b>                   | Anti-ARL2, ADP-ribosylation factor-like protein 2                                                                                                                                                                                                                                                                   |

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**