

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human SCARB2/CD36L2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12045
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human SCARB2/CD36L2 (Gly46-Asp418).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-CD36, 85 kDa lysosomal membrane sialoglycoprotein, CD36 antigen-like 2, CD36L2, SCARB2, Lysosome membrane protein 2, LIMP II, Scavenger receptor class B member 2, Lysosome membrane protein II, LIMPII, LIMP2, LGP85