

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human POU3F2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12051
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human POU3F2 (Thr266-Gln351).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Brain-specific homeobox/POU domain protein 2, POU domain, class 3, transcription factor 2, Brain-2, OTF7, Oct-7, Nervous system-specific octamer-binding transcription factor N-Oct-3, Brn-2, Octamer-binding protein 7, OCT7, Octamer-binding transcripti