

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ZFP36L1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12083
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human ZFP36L1 (Arg121-Arg195).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ERF-1, ZFP36-like 1, TIS11B, BRF1, ZFP36L1, Butyrate response factor 1, Zinc finger protein 36, C3H1 type-like 1, ERF1, RNF162B, BERG36, TPA-induced sequence 11b, mRNA decay activator protein ZFP36L1, EGF-response factor 1