

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Anti-Human MOB1B Polyclonal Antibody |
| PRODUCT REFERENCE | ARO-A12089 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | |
| TARGET | E. coli - derived recombinant Human MOB1B (Ser23-Arg216). |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Oryctolagus cuniculus (Rabbit)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-MOBKL1A, MOB kinase activator 1B, MOB4A, MOB1B, Mob1 homolog 1A, Mob1B, Mob1A, Mps one binder kinase activator-like 1A |