

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human MAD1L1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12097
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human MAD1L1 (His544-Ala718).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-HsMAD1, MAD1-like protein 1, hMAD1, Mitotic arrest deficient 1-like protein 1, Tax-binding protein 181, MAD1, TXBP181, Mitotic checkpoint MAD1 protein homolog, MAD1L1, Mitotic spindle assembly checkpoint protein MAD1